




Quality system certified
to ISO 9001:2008



DENTAL COMPRESSORS

Catalogue 2015



(ENG) NARDI COMPRESSORI was founded in 1965 by the Ghiotto family in Montecchio Maggiore (Vicenza). The symbol  indicates the union of the two names Nardi and Ghiotto owners of the company and leaders in the design and construction of high-tech compressors. In 2008 ongoing investments to improve production processes led to the purchase of a facility of 5000 m² with modern and efficient departments. In 2010, the company was certified UNI EN ISO 9001:2008 with the certifying body TÜV which is among the most qualified in the world. Always careful about production and technology, NARDI COMPRESSORI is proud of having a fully ITALIAN production with testing of each product individually. Today NARDI COMPRESSORI is known in over 100 countries around the world with a

range of compressors from 8 bar up to over 400 bar and a production of over 30,000 units produced each year.

(ITA) NARDI COMPRESSORI è stata fondata nel 1965 dalla famiglia Ghiotto a Montecchio Maggiore (VI). Il simbolo  indica l'unione dei due nomi NARDI e GHIOTTO titolari della azienda, leader nella progettazione e costruzione di compressori ad alto contenuto tecnologico. I continui investimenti atti a migliorare i processi produttivi hanno portato nel 2008 all'acquisto dello stabilimento di 5000 m³ con reparti efficienti e moderni. Nel 2010 l'azienda si è certificata UNI EN ISO 9001:2008 con l'ente certificante TÜV fra i più qualificati al mondo. Sempre attenta alla produzione e alla tecnologia, la NARDI COMPRESSORI vanta di una produzione completamente ITALIANA con collaudi di ogni singolo articolo. Oggi è apprezzata in più di 100 paesi in tutto il mondo con una gamma di compressori che va dai 8 bar fino a oltre i 400 bar e una produzione che supera i 30.000 pezzi prodotti ogni anno.




(ENG) Nardi Compressori S.r.l. has designed and developed a full line of compressors for Dental Industry as a complement to a very wide range of products that covers various market niches. Its internal organization and its dynamic structure allows Nardi Compressori to respond swiftly to market trends and meet customer requirements.

Nardi's Dental compressors are designed and manufactured in-house by our expert technicians with many years of experience in this sector.

Today Nardi Compressori is able to produce oil free Dental compressors which are modern, strong and silent. Dental compressors' tanks are painted inside with epoxy electrostatic powder paint to avoid corrosion.

All compressors are equipped with efficient filtering system that absorbs humidity and dries the air. The filter are made with silica gel that guarantees the air will remain moisture free.



(ITA) La ditta Nardi Compressori S.r.l. ha progettato e sviluppato una linea di compressori dentali a completamento di una vasta gamma di prodotti che coprono diverse nicchie di mercato.

L'organizzazione interna e la struttura dinamica della ditta, le consentono di rispondere velocemente e tempestivamente alle esigenze del mercato ed alle necessità dei clienti.

I compressori dentali sono progettati e prodotti direttamente e totalmente in fabbrica da esperti tecnici che


vantano molti anni di esperienza nel settore.

Oggi Nardi Compressori è in grado di produrre compressori

dentali moderni, resistenti e silenziosi.

I serbatoi dei compressori sono tutti verniciati internamente con polvere epossidica allo scopo di evitarne la corrosione.

Tutti i compressori hanno un efficiente sistema di filtraggio che assorbe l'umidità. Il filtro è costituito da gel di silicio che garantisce all'aria di rimanere secca ed inodore.



Packing

Imballo

Weight with box (ESPRIT 60) : 26 Kg - 57,3 Lbs
 Dimension box : 49 x 49 x 72 (0,017m³)
 Quantity per PALLET (110x110) : 8 pz.




 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
ESPRIT 60/4	1	X	X
ESPRIT 65/4	1	X	X
ESPRIT 60/2	1	X	X
ESPRIT 65/2	1	1	X

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complete	N°	l/min. - C.F.M. - m³/h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
ESPRIT 2V 60/4 25L	25	X	2	105 - 3,7 - 6,3 (50 Hz)	8 - 114	0,75 - 0,5	1400	52	23 - 51	41 x 45 x 52
ESPRIT 2V 65/4 25L	25	X	2	125 - 4,4 - 7,5 (50 Hz)	8 - 114	0,9 - 0,65	1400	53	23,5 - 52	41 x 45 x 52
ESPRIT 2V 60/2 25L	25	X	2	210 - 7,4 - 12,6(50 Hz)	8 - 114	0,75 - 0,5	2800	62	23 - 51	41 x 45 x 52
ESPRIT 2V 60/2 25L	25	✓	2	210 - 7,4 - 12,6(50 Hz)	8 - 114	0,75 - 0,5	2800	62	23 - 51	58,5 x 45 x 52
ESPRIT 2V 65/2 25L	25	X	2	250 - 8,8 - 15 (50 Hz)	8 - 114	0,9 - 0,65	2800	63	23,5 - 52	41 x 45 x 52
ESPRIT 2V 65/2 25L	25	✓	2	250 - 8,8 - 15 (50 Hz)	8 - 114	0,9 - 0,65	2800	63	23,5 - 52	58,5 x 45 x 52

With automatic thermal protection

Con protezione termica automatica

EXTREME 2D 25 Liters



Packing Imballo

Weight with box (EXTREME 2D 1,5 cv) :
44 Kg - 97 Lbs
Dimension box : 51,0 x 56,5 x 71 (0,021m³)
Quantity per (110x110) PALLET : 8 pz.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complite	N°	l/min. - C.F.M. - m³/h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME 2V 25L 0,75 CV	25	X	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	30,5 - 67	47 x 43 x 66,5
EXTREME 2D 25L 0,75 CV	25	✓	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	38,5 - 85	52,5 x 47 x 66,5
EXTREME 2V 25L 1,50 CV	25	X	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	33 - 73	47 x 43 x 66,5
EXTREME 2D 25L 1,50 CV	25	✓	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	41 - 90,5	52,5 x 47 x 66,5
EXTREME 2V 25L 2,00 CV	25	X	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	35,5 - 78	47 x 43 x 66,5
EXTREME 2D 25L 2,00 CV	25	✓	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	44 - 97	52,5 x 47 x 66,5

With automatic thermal protection

Con protezione termica automatica

MADE IN ITALY BY NARDI COMPRESSORI

Page 05



Packing Imballo

Weight with box (EXTREME 2D 1,5 cv) :
51 Kg - 112,4 Lbs
Dimension box : 54 x 57,5 x 90 (0,027m³)
Quantity per (110x110) PALLET : 8 pz.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complite	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME 2V 50L 0,75 CV	50	X	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	36,5 - 81	47 x 43 x 83
EXTREME 2D 50L 0,75 CV	50	✓	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	45 - 99	60 x 47 x 83
EXTREME 2V 50L 1,50 CV	50	X	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	39,5 - 87	47 x 43 x 83
EXTREME 2D 50L 1,50 CV	50	✓	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	48 - 106	60 x 47 x 83
EXTREME 2V 50L 2,00 CV	50	X	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	42,5 - 94	47 x 43 x 83
EXTREME 2V 50L 2,00 CV	50	✓	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	51 - 112	60 x 47 x 83
EXTREME 2V 50L 2,50 CV	50	X	4	480 - 17 - 29 (50 Hz)	8 - 114	2,50 - 1,84	1400	77	52 - 115	64 x 47 x 89
EXTREME 2D 50L 2,50 CV	50	✓	4	480 - 17 - 29 (50 Hz)	8 - 114	2,50 - 1,84	1400	77	60 - 132	64 x 47 x 89

With automatic thermal protection

Con protezione termica automatica

EXTREME 3D 50 Liters



Packing

Imballo

Weight with box (EXTREME 2D 1,5 cv) :

53 Kg - 97 Lbs

Dimension box : 86,5 x 43 x 74 (0,021m³)

Quantity per (110x110) PALLET : 8 pz.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complite	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME 3V 50L 0,75 CV	50	X	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	37 - 82	84 x 31 x 69
EXTREME 3D 50L 0,75 CV	50	✓	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	70	44,5 - 98	84 x 41 x 69
EXTREME 3V 50L 1,50 CV	50	X	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	40 - 88	84 x 31 x 69
EXTREME 3D 50L 1,50 CV	50	✓	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	71	47,5 - 105	84 x 41 x 69
EXTREME 3V 50L 2,00 CV	50	X	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	43 - 95	84 x 31 x 69
EXTREME 3D 50L 2,00 CV	50	✓	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	50,5 - 112	84 x 41 x 69

With automatic thermal protection

Con protezione termica automatica

MADE IN ITALY BY NARDI COMPRESSORI

Page 07



Packing

Imballo

Weight with box (EXTREME SD 1,5 cv) :

78 Kg - 172 Lbs

Dimension box : 73 x 47 x 79 (0,024m³)

Quantity per EUR PALLET : 4 pz.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 0,75	1	X	X
EXTREME 1,50	2	1	1
EXTREME 2,00	2	2	1

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complite	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME SV 30L 0,75 CV	30	X	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	63	65,5 - 145	65 x 44 x 76
EXTREME SD 30L 0,75 CV	30	✓	2	200 - 7,0 - 12 (50 Hz)	8 - 114	0,75 - 0,5	1400	63	72,5 - 159	71 x 44 x 76
EXTREME SV 30L 1,50 CV	30	X	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	63	68 - 150	65 x 44 x 76
EXTREME SD 30L 1,50 CV	30	✓	2	240 - 8,5 - 14,4 (50 Hz)	8 - 114	1,50 - 1,1	1400	63	75 - 165	71 x 44 x 76
EXTREME SV 30L 2,00 CV	30	X	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	64	71,5 - 158	65 x 44 x 76
EXTREME SD 30L 2,00 CV	30	✓	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	64	78,5 - 173	71 x 44 x 76

With automatic thermal protection

Con protezione termica automatica

EXTREME 4D 90 Liters



Packing

Imballo

Weight with box (EXTREME 4D) :

86 Kg - 190 Lbs

Dimension box : 113 x 57 x 84 (0,054m³)

Quantity per EUR PALLET : 2 pz.

 Dental unit	Pneumatic		
	Electric-Pneumatic		
	Electric		
EXTREME 2,00	2	2	1
EXTREME 2,50	5	3	2

Model Compressor	Tank	Dryer	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Essicat.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	Complite	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME 4V 90L 2,00 CV	90	✗	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	71,5 - 158	110 x 54 x 81
EXTREME 4V 90L 2,00 CV	90	✓	2	400 - 14,1 - 24 (50 Hz)	8 - 114	2,00 - 1,5	1400	73	71,5 - 158	110 x 59 x 81
EXTREME 4V 90L 2,50 CV	90	✗	4	480 - 17 - 29 (50 Hz)	8 - 114	2,50 - 1,84	1400	77	71 - 159	110 x 54 x 81
EXTREME 4D 90L 2,50 CV	90	✓	4	480 - 17 - 29 (50 Hz)	8 - 114	2,50 - 1,84	1400	77	80 - 177	110 x 59 x 81


With automatic thermal protection

Con protezione termica automatica


MADE IN ITALY BY NARDI COMPRESSORI

Page 09

*F.A.D




ESPRIT / 60



Model Compressor	Tank	Cylinder	Air Intake	Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Serb.	Cilindri	Aria Aspirata	Pressione	Motore	Giri	Rumore	Peso	Dimensioni
Type	Liters	N°	l/min. - C.F.M. - m³/h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
ESPRIT UNIT 60/4	3	2	105 - 3,7 - 6,3 (50 Hz) 125 - 4,4 - 7,5 (60 Hz)	8 - 114	0,75 - 0,5	1400 1700	48	8 - 17	29 x 22 x 14
ESPRIT UNIT 60/2	3	2	210 - 7,4 - 12,6 (50 Hz) 250 - 8,8 - 15 (60 Hz)	8 - 114	0,75 - 0,5	2800 3300	55	8 - 17	29 x 22 x 14

*F.A.D



ESPRIT / 65





Model Compressor	Cylinder	Air Intake	Max Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Cilindri	Aria Aspirata	Pressione mass.	Motore	Giri	Rumore	Peso	Dimensioni
Type	N°	l/min. - C.F.M. - m³/h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
ESPRIT UNIT 65/4	2	125 - 4,4 - 7,5 (50 Hz) 150 - 5,3 - 9 (60 Hz)	8 - 114	0,90 - 0,65	1400 1700	49	8,5 - 19	29 x 22 x 14
ESPRIT UNIT 65/2	2	250 - 8,8 - 15 (50 Hz) 300 - 10,5 - 18 (60 Hz)	8 - 114	0,90 - 0,65	2800 3300	56	8,5 - 19	29 x 22 x 14



With automatic thermal protection - * F.A.D. values have a tolerance $\pm 5\%$

Con protezione termica automatica - * I valori indicati hanno una tolleranza $\pm 5\%$

PUMP UNIT



Model Compressor	Cylinder	Air Intake	Max Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Cilindri	Aria Aspirata	Pressione mass.	Motore	Giri	Rumore	Peso	Dimensioni
Type	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME UNIT 0,75 CV	2	200 – 7,0 – 12 (50 Hz) 240 – 8,5 – 14,4 (60 Hz)	10 – 145	0,75 – 0,55	1400 1700	72	14 – 31	39 x 28 x 26
EXTREME UNIT 1,50 CV	2	240 – 8,5 – 14,4 (50 Hz) 290 – 10,1 – 17,2 (60 Hz)	10 – 145	1,50 – 1,1	1400 1700	72	16,5 – 36	39 x 28 x 26




Model Compressor	Cylinder	Air Intake	Max Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Cilindri	Aria Aspirata	Pressione mass.	Motore	Giri	Rumore	Peso	Dimensioni
Type	N°	l/min. - C.F.M. - m ³ /h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME UNIT 2 CV 1400	2	350 – 12,4 – 21 (50 Hz) 420 – 14,8 – 25 (60 Hz)	10 – 145	2,00 – 1,5	1400 1700	73	20 – 44	43 x 28 x 29

With automatic thermal protection - * F.A.D. values have a tolerance $\pm 5\%$

Con protezione termica automatica - * I valori indicati hanno una tolleranza $\pm 5\%$

MADE IN ITALY BY NARDI COMPRESSORI

Page 11



EXTREME 70



Model Compressor	Cylinder	Air Intake	Max Pressure	Power	Revolution	Noise	Weight	Dimension
Modello Compressore	Cilindri	Aria Aspirata	Pressione mass.	Motore	Giri	Rumore	Peso	Dimensioni
Type	N°	l/min. - C.F.M. - m³/h	Bar - Psi	HP - kW	Rpm	dB	Kg - Lbs	cm (LxWxH)
EXTREME UNIT 2,50 CV	4	480 - 17 - 29 (50 Hz) 580 - 20 - 35 (60 Hz)	10 - 145	2,50 - 1,84	1400 1700	73	29 - 64	61 x 28 x 31

DRYER 400



Model	Air filtration	Max Pressure	Weight	Dimension
Modello	Aria filtrata	Pressione mass.	Peso	Dimensioni
Type	l/min. - C.F.M. - m³/h	Bar - Psi	Kg - Lbs	cm (LxWxH)
DRYER 400	400 - 14 - 24	10 - 145	7 - 15,5	13 x 26 x 51

With automatic thermal protection - * F.A.D. values have a tolerance $\pm 5\%$

Con protezione termica automatica - * I valori indicati hanno una tolleranza $\pm 5\%$

General terms



General terms of sale (ENGLISH)

Standard Delivery Fornitura standard

Standard Volt/Hertz/Phase : 230/50/1
Other : 100⇒240 Volt/50-60 Hz/1 Phase

STANDARD COLOR : WHITE
Other colors min. delivery : 50 pcs

Standard plug :
SCHUKO

Standard fitting :
1/4" G BALL VALVE

- This catalogue cancels and substitutes the previous versions.
- Once order is confirmed all terms of sale are implicitly accepted.
- The goods are always delivered ex works and at the buyer's risk.
- Any possible complaints about the packing or the delivery are accepted only if notified within 8 days from the date of receipt.
- Our production is warranted for 24 months from the delivery of the goods. Warranty is excluded for improper use and consequent damages. Warranty is guaranteed only for manufacturing fault.
- Transport costs are always and completely at customer charge.
- Delivery terms are always to be considered as approximate. Shortage of all or part of the goods on delivery can be claimed only if accompanied by written terms agreed upon. Otherwise redraft cannot be claimed.
- All data and technical aspects cited in this catalogue are not binding and may be altered without notice due to technical or manufacturing necessities.
- The invoice referring to the delivered goods is to be paid according to the terms and in the place agreed. Any delay in paying it will enforce late charges.
- In case of disputes, jurisdiction will be the court of Vicenza – Italy.
- This catalogue belongs exclusively to the firm NARDI COMPRESSORI situated in Montecchio Maggiore (VI) Italy; reproduction or copying of the goods without the firm's permission is against the law.

Condizioni generali di vendita (ITALIANO)

- Il presente catalogo annulla e sostituisce il precedente.
- Con il conferimento dell'ordine si considerano accettate tutte le condizioni di vendita.
- La merce si considera franco nostra fabbrica. La merce viaggia sempre a rischio e pericolo dell'acquirente che ha l'obbligo di ritirarla.
- Non si accettano reclami su imballi o spedizioni se non fatti entro 8 giorni dalla consegna della merce.
- I nostri prodotti sono coperti di garanzia per 24 mesi dalla consegna della merce. Non sono garantiti i danni provocati da usura, cattivo uso e non ci assumiamo nessuna responsabilità per danni diretti o indiretti derivanti da tale cattivo uso. La garanzia è valida solo per difetti di fabbricazione.
- Tutti i costi di trasporto sono sempre e completamente a carico del cliente.
- I termini di consegna non sono mai impegnati. La mancata consegna di parte o tutta l'ordinazione non può dar luogo ad alcuna rivalsa di danni verso la nostra ditta, salvo diverse attuazioni scritte e regolarmente accettate.
- Tutti i dati e le indicazioni tecniche citate nel presente catalogo non sono impegnativi e potranno essere modificati in caso di necessità tecniche o produttive, senza alcun preavviso.
- Per pagamenti effettuati oltre la scadenza pattuita saranno addebitati l'interessi moratori al tasso corrente bancario.
- In caso di controversia è riconosciuta la sola competenza del tribunale di Vicenza.
- Questo catalogo è di proprietà esclusiva della ditta Nardi Compressori Srl con sede a Montecchio Maggiore (VI) e senza il consenso della stessa ditta non può essere né copiato, né riprodotto anche solo parzialmente sotto pena delle sanzioni di legge in materia.


This catalogue replaced the previous one

Questo catalogo sostituisce i precedenti

MADE IN ITALY BY NARDI COMPRESSORI

Page 13

Nardi Compressori ITALY



Connection



Railway station **VICENZA**

15 Km



Airport **VERONA** "Villafranca"

VRN

55 Km



Airport **VENICE** "Marco Polo"

VCE

83 Km



Airport **MILAN** "Orio al Serio"

BGY

150 Km



Airport **BOLOGNA** "Marconi"

BGY

170 Km



Airport **MILAN** "Linate"

LIN

200 Km




Airport **MILAN** "Malpensa"

MLP

230 Km

PRODUCTION



A complete line of oil free low pressure compressors with maximum working pressure of 10 Bar plus a new line of 30 Bar compressors. These compressors have different capacity vessels up to a maximum of 90 liters. Power includes single phase and 3-phase AC electric motors, DC electric motors (battery) as well as internal combustion engines.



High pressure compressors up to 450 Bar for the filling of cylinders for breathing air, used often for scuba diving and by fire stations for fire and rescue operations as well as for military applications.



High pressure compressors and boosters up to 450 Bar. Technical gas compressors with high degree of purity for industrial use, the range is completed by the supply of high pressure cylinders for storage in different capacity and configuration.



High pressure compressors for refueling of vehicles that run on natural gas. Compressors in this line have been designed for industrial and domestic use, creating simple and economic small refueling stations for private use using gas supplied through the network lines.





Nardi Compressori Srl
Via Marco Polo, 2 36075
Montecchio Maggiore (VICENZA) ITALY
Tel. +39.0444.159111 - 497514 -694885
Fax. +39.0444.159122
Mail. info@nardicompressori.com



www.nardicompressori.com